

WATER BATH 6L

WATER BATHS

MAIN FEATURES:

- > Small Bench Top Water Bath
- > Stainless steel chamber and cover
- > Powerful heater controlled by LED digital controller
- > Precise temperature control by hydraulic thermostat

Digital Water Bath is solidly constructed with a stainless steel, seamless chamber that resists corrosion. A stainless gable cover is included to minimize evaporation. The optional non-toxic (non-mercury) thermometer clips onto the gable cover for accurate temperature monitoring. The chamber can accommodate plastic tube racks.



TECHNICAL SPECIFICATIONS

Model	Digital Water Bath 6l
Voltage	230 V / 50-60 Hz
Power	500VA
Temperature range	Ambient + 5°C to 99°C
Controls	PID Control
Temperature accuracy at 37°C	±0,5°C
Temperature stability at 37°C	±0,3°C
Resolution	0,1°C
Chamber volume	6 liters
Chamber maximum load	9 kg
Net weight	8,6 kg
Dimensions	380 x 290 x 330 mm
Chamber Dimensions	300 x 150 x 150 mm

ORDERING INFORMATION

WATER BATH 6L, 230V / 50-60Hz	785.500
-------------------------------	---------

We reserve the right to alter specification details etc. without prior notice or liability!

V:11-19

DOMEL

DOMEL, d.o.o., Otoki 21, 4228 Železniki
BU Laboratory Systems, Na Plavžu 79, 4228 Železniki

Phone: +386 4 5117 525
E-mail: info.labs@domel.com
WEB: <http://labs.domel.com>

