

ACCESSORIES EMBEDDING



SD 2800 Slide Dryer

IVD CE



The unit

The SD2800 is a heated plate designed to dry histological slides.

The Teflon coating allows easy cleaning.

The SD 2800 is provided with a control panel to set and monitor the temperature of the heating surface.

Is equipped with overheating protection.

The area temperature is maintained by a microprocessor.

Technical features:

Temperature range

Heated plate: ambient ~ 75°C increments 1°C

Dimensions and weight

Width: 324 mm

Depth: 272 mm

Height: 83 mm

Weight: 2,75 kg

Electrical data

Power voltages: 220V - 50 Hz;

Maximum power draw: 300 V/A



Product	Size	Packaging	Code
Slide Dryer	Unit	1 pc	SD2800

ACCESSORIES EMBEDDING



Differences between instruments

Through a wide range of accessories for the Histo-Line embedding it is possible to equip the pathology laboratories in different ways depending on the space available and the operators' requests.



Embedding Center Accessories

Peculiarity	TEC2500	WB2900	WB2800	SD2800
Adjustable temperature	X	X	X	X
Programmable	X			
Teflon surface	X	X	X	X
Microprocessor control	X	X	X	X
Parts removable	X	X		